

## mikrozyd® wipes

Safety on all surfaces.



### mikrozyd® AF wipes

The classic choice among rapid disinfection preparations. The alcohol-based wipes are effective within 1 minute against most relevant microbes. The alcohol film dries very quickly and leaves no residues behind. Material compatibility is good, however stress cracks may form in plastics.

**Shelf life\*:** 3 months | **Active ingredients:** Ethanol, Propan-1-ol

Article	Article number	Ø Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	109203	Ø 14 x 18 cm	150	24 g/m <sup>2</sup>	0.3 m <sup>2</sup>	
Tub Refill	70002062	Ø 14 x 18 cm	150	24 g/m <sup>2</sup>	0.3 m <sup>2</sup>	
Jumbo Tub	70003093	Ø 20 x 20 cm	220	24 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Jumbo Tub Refill	70003094	Ø 20 x 20 cm	220	24 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Softpack	70000808	Ø 20 x 20 cm	50	50 g/m <sup>2</sup>	0.5 m <sup>2</sup>	

### mikrozyd® sensitive wipes

The specialist for sensitive surfaces. There are many sensitive materials in a hospital and for this reason, mikrozyd® sensitive wipes combine a high degree of efficacy with excellent material compatibility. Even ultrasound heads\*\* can be disinfected with mikrozyd® sensitive wipes.

**Shelf life\*:** 3 months | **Active ingredient:** QAC perfume-free

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	70003103	Ø 14 x 18 cm	120	20 g/m <sup>2</sup>	0.3 m <sup>2</sup>	
Jumbo Tub	70003104	Ø 20 x 20 cm	200	20 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Jumbo Tub Refill	70003109	Ø 20 x 20 cm	200	20 g/m <sup>2</sup>	0.5 m <sup>2</sup>	
Softpack	70003101	Ø 20 x 20 cm	50	50 g/m <sup>2</sup>	1 m <sup>2</sup>	
Softpack	70003108	Ø 20 x 20 cm	100	50 g/m <sup>2</sup>	1 m <sup>2</sup>	
Maxi-Softpack	70003102	Ø 25 x 25 cm	80	50g/m <sup>2</sup>	1.5m <sup>2</sup>	

### mikrozyd® universal wipes

The all-round wipes – also for disinfecting sensitive communication devices. Smartphones and tablets have nearly constant contact with hands. They can be safely treated with mikrozyd® universal wipes without damage to the materials.

**Shelf life\*:** 1 month | **Active ingredients:** Ethanol, Propan-2-ol perfume-free

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Softpack	70000038	Ø 20 x 20 cm	100	50 g/m <sup>2</sup>	1 m <sup>2</sup>	
Maxi-Softpack	70001061	Ø 25 x 25 cm	80	50g/m <sup>2</sup>	1.5m <sup>2</sup>	

### mikrozyd® PAA wipes

The powerhouse for cases of outbreak. With the unique oxidising power of peracetic acid, mikrozyd® PAA wipes have a sporicidal effect and are also effective against C. difficile in just 5 minutes.

**Shelf life\*:** 28 days | **Active ingredient:** Peracetic acid perfume-free

Article	Article number	Wipe size	Wipe pieces	Grammage	Surface coverage Ø	Efficacy
Tub	70003205	Ø 20 x 30 cm	50	40 g/m <sup>2</sup>	3.5 m <sup>2</sup>	

Use disinfectants safely. Always read the label and product information before use.

\* The shelf life after the container is opened is marked on articles with the following symbol:

\*\* Not suitable for final disinfection of semicritical medical devices.